



MIDDLE EAST

FOR NON-WOVEN CLOTHES

Surgery Packs

Book Now

sales@middleeastmed.com

WWW.MIDDLEEASTMED.COM

Materials used chart

SMS (Spunbond-Meltblown-Spunbond) :

- عبارة عن ثلاث طبقات من ألياف البولي بروبيلين (PP)، توفر قوة جيدة مع مقاومة عالية للسوائل.
- يُستخدم عادةً في العباءات الجراحية، الأغذية الطبية، والأقنعة.
- Consists of three layers of polypropylene (PP) fibers, providing good strength with high liquid resistance.
- Used in surgical gowns, medical covers, and masks.

SMS REINFORCED :

- عبارة عن SMS لكن عليه طبقة إضافية أو منطقة محددة لدعم الحماية (خصوصاً في مناطق الصدر والأكمال).
- يعزز مقاومة الاختراق من السوائل ويزيد الحماية في المناطق الأكثر عرضة للتلوث.
- It is SMS with an additional layer or a specific area to enhance protection (especially in the chest and sleeve areas).
- Enhances liquid penetration resistance and increases protection in areas more prone to contamination.

PE (Polyethylene) :

- طبقة بلاستيكية غير منسوجة، مقاومة ممتازة للسوائل.
- تُستخدم عادةً في الأغذية والعباءات التي تتطلب حاجزاً عالياً ضد السوائل.
- A nonwoven plastic layer with excellent liquid resistance.
- Commonly used in covers and gowns that require a high barrier against liquids.

spunlace non woven fabric :




- قماش غير منسوج مصنوع عن طريق تشابك الألياف بالماء عالي الضغط (hydroentanglement).
- ناعم الملمس وبشبه القماش المنسوج، يُستخدم في الشراشف والملابس التي تتطلب راحة عالية.
- Nonwoven fabric made by entangling fibers with high-pressure water (hydroentanglement).
- Soft to the touch and similar to woven fabric, used in sheets and garments requiring high comfort.

Adhesive :

- يشير إلى مناطق بها طبقة لاصقة تُستخدم لتثبيت الأقمشة أو تأمين مواضع معينة (مثل الشرائط الجراحية).
- Refers to areas with an adhesive layer used to secure fabrics or fix specific positions (such as surgical tapes).

Opsite :


- فيلم شفاف لاصق مقاوم للماء وغالباً يستخدم لتغطية الجروح أو مناطق صغيرة لحمايتها.
- يسمح بالرؤية من خلاله ويحافظ على الحاجز ضد السوائل والبكتيريا.
- Transparent adhesive waterproof film, commonly used to cover wounds or small areas for protection.
- Allows visibility through it while maintaining a barrier against fluids and bacteria.






Contents of the Surgical pack A:

Code	Name	Size	Q
ME-S122 A	Surgical Drape A	300×160 cm	1
ME-S126 C1	Adhesive Drape C	100×80 cm	4
ME-S125	Mayo Stand Cover	140×80 cm	1
ME-S124 A2	Instrument Trolley Cover A2	180×150 cm	1



Contents of the Surgical pack A1:

Code	Name	Size	Q
ME-S122 A	Surgical Drape A	300×160 cm	1
ME-S126 C1	Adhesive Drape C1	100×80 cm	4
ME-S125	Mayo Stand Cover	140×80 cm	1
ME-S124 A2	Instrument Trolley Cover A2	180×150 cm	1
ME-S216 D	Hand Towel D	20×30cm	4
ME-S221	Diathermy Cover	30×40cm	1




MIDDLE EAST
FOR NON-WOVEN CLOTHES






Contents of the Surgical pack A2:

Code	Name	Size	Q
ME-S122 A	Surgical Drape A	300×160 cm	1
ME-S126 C1	Adhesive Drape C1	100×80 cm	4
ME-S125 A	Mayo Stand Cover A	140×80 cm	1
ME-S124 A2	Instrument Trolley Cover A2	180×150 cm	1
ME-S127 A4	Utility Drape A4	80×80 cm	2
ME-S221	Diathermy Cover	30×40 cm	1



Contents of the Surgical pack A3:

Code	Name	Size	Q
ME-S122 A2	Surgical Drape A2	200×160 cm	1
ME-S126 C	Adhesive Drape C	80×80 cm	6
ME-S217	Mayo Cover SMS	75×120 cm	1



MIDDLE EAST
FOR NON-WOVEN CLOTHES






Contents of the Surgical pack S:


Code	Name	Size	Q
ME-S122 A	Surgical Drape A	300×160 cm	1
ME-S126 C1	Adhesive Drape C1	100×80 cm	4
ME-S101 S02	Surgical Gown	XXL	3
ME-S125	Mayo Stand Cover	140×80 cm	1
ME-S124 A5	Instrument Trolley Cover A5	180×150 cm	1
ME-S216 D	Hand Towel D	20 x 30 cm	3



ME-S126 C1




ME-S122 A



ME-S124 A5



ME-S125



ME-S101 S02




ME-S216 D

Contents of the Surgical pack S1:


Code	Name	Size	Q
ME-S122 A	Surgical Drape A	300×160 cm	1
ME-S126 C1	Adhesive Drape C1	100×80 cm	4
ME-S101 S02	Surgical Gown	XXL	2
ME-S125	Mayo Stand Cover	140×80 cm	1
ME-S124 A5	Instrument Trolley Cover A5	180×150 cm	1
ME-S216 D	Hand Towel D	20 x 30 cm	3



ME-S126 C1



ME-S122 A



ME-S124 A5



ME-S125



ME-S101 S02




ME-S216 D






Contents of the Surgical pack B:


Code	Name	Size	Q
ME-S122 A1	Surgical Drape A1	250X160 cm	1
ME-126 C1	Adhesive drape	100X80 cm	4
ME-125	Mayo Stand Cover	140X80 cm	1
ME-124 A5	Instrument Trolley Cover	160X150 cm	1



ME-S126 C1



ME-S122 A1



ME-S124 A5



ME-S125

Contents of the Surgical pack C1:

Code	Name	Size	Q
ME-122 A2	Surgical Drape A2	200×160 cm	1
ME-S101 A04	Surgical Gown	Large	2
ME-S126 C1	Adhesive Drape C1	80×80 cm	4
ME-S127 A4	Utility Drape A4	80×80 cm	1



ME-S126 C1



ME-S122 A2



ME-S101 A04




ME-S127 A4

Contents of the Surgical pack C2:

Code	Name	Size	Q
ME-S122 A2	Surgical Drape A2	200×160 cm	1
ME-S101 A04	Surgical Gown	Large	3
ME-S126 C	Adhesive Drape C	80×80 cm	4
ME-S127 A4	Utility Drape A4	80×80 cm	1



ME-S126 C



ME-S122 A2



ME-S101 S02
ME-S101 A04



ME-S127 A4






Contents of the Surgical pack C3:


Code	Name	Size	Q
ME-S122 A2	Surgical Drape A2	200×160 cm	1
ME-S101 A04	Surgical Gown	Large	3
ME-S126 C	Adhesive Drape C	80×80 cm	4
ME-S125 A1	Mayo Stand Cover A 1	100×80 cm	1
ME-S127 A4	Utility Drape A4	80×80 cm	1



ME-S126 C



ME-S122 A2



ME-S101 A04



ME-S125 A1




ME-S127 A4

Contents of the Surgical pack C4:

Code	Name	Size	Q
ME-S122 A2	Surgical Drape A2	200×160 cm	1
ME-S101 A04	Surgical Gown	Large	2
ME-S126 C	Adhesive Drape C	80×80 cm	4
ME-S125 A1	Mayo Stand Cover A 1	100×80 cm	1
ME-S127 A4	Utility Drape A4	80×80 cm	1



ME-S126 C



ME-S122 A2



ME-S101 A04



ME-S125 A1




ME-S127 A4






Contents of the Universal Pack A:


Code	Name	Size	Q
ME-125	Mayo Stand Cover	140×80 cm	1
ME-124 A2	Instrument Trolley Cover A2	150×180 cm	1
ME-126 A1	Adhesive Drape A1 (Leg)	180×160 cm	1
ME-S102 S02	Full Reinforced Gown	XXL	2
ME-126 A3	Adhesive Drape A3 (Head)	250×160 cm	1
ME-126 C1	Adhesive Drape C1	100×80 cm	1
ME-216 A	Hand Towel A	30×40 cm	2
ME-220 A	OP-Tape	10×50 cm	2



ME-126 A3



ME-126 A1




ME-124 A2




ME-125



ME-S102 S02



ME-126 C1




ME-220 A



ME-216 A

Contents of the Universal Pack B:

Code	Name	Size	Q
ME-125	Mayo Stand Cover	140×80 cm	1
ME-124 A2	Instrument Trolley Cover A2	150×180 cm	1
ME-126 A1	Adhesive Drape A1 (Leg)	180×160 cm	1
ME-S101 S02	Surgical Gown	XXL	2
ME-126 A3	Adhesive Drape A3 (Head)	250×160 cm	1
ME-126 C1	Adhesive Drape C1	100×80 cm	1
ME-216 A	Hand Towel A	30×40 cm	2
ME-220 A	OP-Tape	10×50 cm	2




ME-126 A3




ME-126 A1




ME-125




ME-124 A2



ME-S101 S02



ME-126 C1



ME-220 A



ME-216 A





MIDDLE EAST
FOR NON-WOVEN CLOTHES

الشرق الاوسط للملابس الغير منسوجه

Beside you at every step so you can provide the best care

For more information or to order, please contact your MIDDLE EAST FOR NON WOVEN sales representative, call +20 11 29192312 or send mail to sales@middleeastmed.com

Explore more Middle East For Non-Woven Clothes products



www.middleeastmed.com



+20 11 29192312